

PANASIL® PUTTY FAST, PUTTY SOFT AND PUTTY



Panasil® Putty is a classic kneadable putty based on A silicone for precision impressions. After more than 40 years of continuous development, the kneading putty is still a reliable quality material.



Easy processing

- Accustomed handling without additional equipment
- Smooth kneading without sticking

Precise and cost effective

- Hardness and elasticity in harmony for easy removal from the mouth
- Extra high dynamic pressure for optimal flow
- Dimensionally accurate
- Outstanding value for money

Always the right product

- Panasil® Putty Soft – reduced final hardness (Shore A 60) for even easier removal from the mouth
- Panasil® Putty Fast – fast setting (end of setting 3:30 min.)
- Panasil® Putty – final hardness Shore A 66 (end of setting 4:00 min.)



Panasil® Putty Fast short intraoral setting time

- + Two-step impressions
- + Foil impressions


Panasil® Putty Soft reduced final hardness




- + One-step putty-wash impressions
- + Functional margin contouring
- + Two-step impressions
- + Foil impressions

Panasil® Putty high dynamic pressure

- + Two-step impressions
- + Foil impressions
- + One-step putty-wash impressions
- + Functional margin contouring

- + Highly recommended
- + Recommended



Panasil®	900 mL Normal pack 1:1 
	1 x 450 mL catalyst paste, 1 x 450 mL base paste, 2 dispensing scoops
Putty Fast 	REF 11141
Putty Soft 	REF 11121
Putty	REF 11101

Follow us:



For more information contact Kettenbach direct **877-532-2123**
website: www.kettenbach-dental.us
website: www.kettenbachusa.com