

**1 IN 4 PATIENTS...**  
**are at risk of having obstructive sleep apnea.**



**The average dental office has 1200 PATIENTS.**  
**A QUARTER OF THESE have sleep apnea!**

## **i** AIRWAY METRICS

### **What is Airway Metrics?**


A system used to obtain an optimum treatment position for oral appliance therapy. Airway metrics helps the clinician identify the optimum protrusive and vertical placement of the mandible which will open the airway in the management of sleep disorders.

### **Why Incorporate this into your practice?**

Your chance to inform patients of a better quality of life and or potentially save their life.  
Strong revenue source \$2,500 to \$5,000 per Mandibular Advancement Device (MAD) dependent of offered therapy.

### **How does Airway Metrics System work?**

Airway Metrics is a system that helps quickly & easily determine the correct treatment position of the vertical and protrusive settings for the bite. Using the Snore Screener and Mandibular Positioning Simulators (MPS) along with your patient's subjective feedback you can identify the correct bite position for your patients optimum airway. Once the correct setting is identified, the recordings are then sent to the lab for an appliance to be fabricated.



The National Institute of Health estimates that between **50 and 70 million Americans** suffer from **chronic sleep disorders**, and that **80 to 90%** of them remain undiagnosed.\*

\*(Chest; 2006; 130; 780-786)



## BENEFITS TO DENTIST

- Designed to promote optimum interaction and patient comfort during the sleep appliance screening process.
- Snore Screener identifies patient's optimum airway zone and determines a comfortable starting position.
- Mandibular Positioning Simulators (MPS) fine tunes the mandibular protrusive and vertical measurements to determine ultimate treatment position.
- Vertical positions of 4, 6, 8, 10 and 12mm allow for maximum adaptability for each patient.
- Incisal position groove for easy positioning.
- Easily read measurement guide.
- Available in 3 colors (blue, amber, black)
- Reusable with cold sterilization.

## Other Kettenbach Products commonly used with Airway Metrics system.

### Impression Material:

- Panasil® Putty Soft & Panasil® initial contact X-Light
- Or Panasil® Tray Soft and Panasil® initial contact X-Light
- Or Silginat® (PVS Alginate Alternative)



### Bite Registration:

- Futar® D Fast



**BUY SMART.  
BUY DIRECT.**

**Subscribe to Kettenbach Corporate YouTube channel for more information on Airway Metrics.**

<https://www.youtube.com/watch?v=Q09FfVOEV-o&t=45s>



JOIN US!

For more information contact Kettenbach direct 877-532-2123  
[www.kettenbachusa.com](http://www.kettenbachusa.com)